



PCR of pieces for ~~...~~ His... | EGR-Cys (PAX 38 → end)

Monday 2/17/03

## FORWARD PRIMER

10 μM

BVPAXFOR (029473-1)

02/13/2003

36 mer

5'-ATA GAA GGA AGA TGC CAC ACA TAC CAG  
GAG ATT GCC1 OD = 29.9 μg = 2701 pmol  
MW = 11084.33 E260 = 370.2 L/mmol-cm  
Tm = 70.1 °C ([Oligo]=250 pM; [Salt]=50 mM)

Yield = 3.6 OD; Purified (RPC) 9723.6 pmol

AMITOF (617)782-9242 • 1-(800)998-4863 • (617) 782-9352(FAX)

## BACKWARD PRIMER

20 μM

BVPAXBAC (029473-2)

02/13/2003

24 mer

5'-CTA GCA GAA GAG CTT GAG GAA GCA

1 OD = 29.9 μg = 4024 pmol  
MW = 7435.95 E260 = 248.5 L/mmol-cm  
Tm = 56.3 °C ([Oligo]=250 pM; [Salt]=50 mM)


Yield = 4.5 OD; Purified (RPC) 18108.0 pmol

AMITOF (617)782-9242 • 1-(800)998-4863 • (617) 782-9352(FAX)

	# 1	# 2	# 3
1	10 × PCR buffer	5 mL	5 mL
2	3 mM dNTPs	5 mL	2 mL
3	50 mM MgCl <sub>2</sub>	1.5 mL	1.5 mL (2 μM)
4	Forward primer (1 μM)	5 mL	5 mL
5	Back primer (1 μM)	2.5 mL	2.5 mL
6	DNA Template (27.5 ng)	0.5 mL	0.5 mL
7	Taq Polymerase* (1 unit/mL)	1 mL	1 mL
8	sterile water (μL)	29.5	32.5
	50 mL	50 mL	50 mL

\* Note: Do not use the platinum Taq - it has proofreaders that will get rid of the A overhangs needed for insertion into the vector. (perhaps just Hi Fi that has proofing)

## PCR products



lane 1: 500 bp molecular ruler

lane 2: PCR conditions 1

lane 3: PCR conditions 2

lane 4: PCR conditions 3

Yea molecular biology!

1 2 3 4

gel: 1.5% agarose in TAE buffer

500 bp molecular ruler (500, 1000, 1500, 2000, etc.)

TOPO Cloning and Transformation of PAX (Cys, 38 → end)

Tuesday 2/18/03

Cloning Reaction (w/ Chemically Competent E. coli)

Fresh PCR Product	2 $\mu$ L (from PCR #2 p. 3-73)
salt soln. (1.2 M NaCl, 0.06 M MgCl <sub>2</sub> ) (from TOPO Kit)	1 $\mu$ L
sterile water	2 $\mu$ L (to 5 $\mu$ L total so far)
TOPO <sup>®</sup> vector	1 $\mu$ L

- mixed gently and incubated 5 min. at room temp. (Flicked w/ finger for mixing, then tapped back down.)
- Reaction was placed on ice until the Transformation.

Transformation Reaction (Following directions for the One Shot<sup>®</sup> TOPO<sup>®</sup> Transformation Rxn. p. 12 TOPO instructions)

Following the transformation 2 LB plates were spread, one with 10  $\mu$ L and one with 50  $\mu$ L from the transformation.


To spread, a spreader sterilized by burning EtOH was used on an antibiotic-containing LB plate. 50  $\mu$ L of SOC were added to each plate before spreading. The plates were incubated overnight at 37° C.

Wednesday 2/19/03


- 5 colonies from each plate were selected for analysis. With a small (pink size) pipet tip, the colony was touched. The tip was then released into a clear snap-like loose top tube containing 3 mL of LB media w/ an antibiotic. The tubes were numbered #1 - #10 and incubated overnight shaking at 37° C.

## Mini Prep and Restriction Digest

Thursday, 2/21/03



- Minipreps on the 10 colonies (#1-10 grown up on the previous page) were performed following the miniprep kit conditions.
- The products (#1-10) of the minipreps were digested with EcoRI. This is an initial screen to determine which colonies have the TOPO vector (either with or w/out the PAX insert) and which colonies have (Debbie's) DNA template used for the PCR rxn.



- The first digest (EcoRI) is simply to distinguish whether or not the colonies contain the DNA template or the TOPO vector in some form. Once this is determined, the colonies containing the TOPO vector will be digested with BamHI. This will determine whether the vector contains the PAX insert and if that insert was ligated in forwards or backwards.