


Repeat of Transformation (w/ new TOP 10F' cells) of → Tuesday, March 4, 2003

Purified PCR product w/ 3'-A overhang addition

- The PCR product cloned into the TOPO vector was transformed into the One shot[®] E. coli cells following the TOPO vector instructions. (Saturday, 3/1/03)
- The transformation product was plated in two different amounts (10 μ L and 50 μ L) onto LB plates w/ carbanecillin. 50 μ L of SOC media was added to aid in spreading.
- (Sunday 3/2/03) A decent number of colonies had grown on the 50 μ L plate and a few on the 10 μ L plate. 10 colonies (#1-10) were selected and grown up in 3mL of LB media w/ carbenicillin overnight, shaking at 37° C.
- ~~10~~ minipreps were done for #1-6 (ran out of ~~spin columns~~ cent. matrixes from eppendorf, so only did 6 (not 10) preps.
- EcoRI digest of miniprep # 1-6 (done as on p. 3-78)



- 1 500 bp marker
- 2 EcoRI digest #4
- 3 #4
- 4 EcoRI digest #3
- 5 #3
- 6 EcoRI digest #2
- 7 #2
- 8 EcoRI digest #1
- 9 #1
- 10 500 bp ladder

- #2 - #6 look like the TOPO vector closed in on itself (2,872 bp)

- #1 looks promising...
(vector + insert should be 4,450 bp. Will do a BamHI digest ~ see p. 3-75 for drawing)

1 2 3 4 5 6 7 8 9 10

BamHI digest (continued from p. 82)

Wednesday, March 5, 2003

BamHI digestion (on # 1-6)


DNA from miniprep
 BamHI 10x rest. buffer
 water
 200 µg/mL BSA soln.
 BamHI enzyme (added
 last and always kept on ice)

8 µL
 2 µL
 9 µL
 1 µL
 0.25 µL

 20 µL

stock (for immediate use)
 30 µL
 15 µL
 135 µL
 3.75 µL

added 12.75 µL to each DNA sample



- In 1 500 bp ladder
- 2 # 1
- 3 BamHI digestion # 1
- 4 # 2
- 5 BamHI digestion # 2
- ← 6 # 3
- 7 BamHI digestion # 3
- 8 # 4
- 9 BamHI digestion # 4
- 10 500 bp ladder

1 2 3 4 5 6 7 8 9 10

- In 1 500 bp ladder
 - * 2 # 1
 - * 3 EcoRI digestion # 1
 - * 4 BamHI digestion # 1
 - 5 # 5
 - 6 BamHI digestion # 5
 - 7 # 6
 - 8 BamHI digestion # 6
 - 9 purified PCR product (p. 3-81)
 - 10 purified PCR w/ A-overhangs added (p. 3-80).
- } # 1 sent for sequencing



1 2 3 4 5 6 7 8 9 10