


Transformation into BL(21)-codon Plus (RP) cells

- Same procedure as for transformation into the RIL cells, described on p. 3-109.
- 100 μ L and 150 μ L of the transformed cells were plated onto LB agar plates with 50 μ g/mL carbenecillin and 50 μ g/mL chloramphenicol (see p. 106 for plate recipe)

RP Pilot Expression:




1 2 3 4 5 6 7 8 9 10

in	1	protein ladder
	2	colony 1 uninduced (w/ 1 M IPTG)
	3	colony 1 induced
	4	" 2 uninduced
	5	" 2 induced
	6	" 3 uninduced
	7	" 3 induced
	8	" 4 uninduced
	9	" 4 induced
	10	" 5 induced

RP Cells Continued:

- 200 mL prep, followed by Ni purification (These were grown from frozen culture overnight in 5 mL. 5 mL were added to 1 L, grown at 37°C to 0.6 OD₆₀₀ and induced w/ 1mM IPTG for 4 hr.



- each lane loaded w/ 20 μ L of a mixture of sample (from 1.5 mL fractions off the Ni column) using 5X SDS loading buffer.

(from colony 3 p. 3-III)

- Prep 2: 1 L grown up to an OD₆₀₀ = 0.5 at 37°C, then placed in a water bath shaker at 16°C for 15 min and induced w/ 1mM IPTG. The 1 L was shaken at 16°C overnight before the cells were spun down and frozen at -80°C.
- The cells were lysed and purified (all 1 L) on a Ni column.

