

Yet another Factor Xa test cleavage (Novabiochem Fact. Xa)

To each of 5 50- μ L eppendorf tubes was added:

- { 6 μ L Factor Xa buffer 10X
- { 24 μ L (\approx 12 μ g) FLAG and GST-purified protein
- { 29 μ L DI water
- 1 μ L Factor Xa (added after 10 μ L of the above mix was removed for gel analysis)

- After the 10X buffer, protein solution and water were added, the tubes were flicked to mix and 10 μ L were added to 2.5 μ L of SDS loading buffer for subsequent gel analysis (as the $t=0$ timepoint).
- Dilutions of the 1 μ g/ μ L neb Factor Xa solution ~~was~~ made to give a 0.5, 0.1, 0.05, 0.025 and 0.0125 dilutions were made. 1 μ L of Fact. Xa was added to each tube as follows:

tube a:	0.2 μ g	d:	0.025 μ g
b:	0.1 μ g	e:	0.0125 μ g
c:	0.05 μ g	f:	0 μ g

- 10- μ L aliquots were removed after 1, 3, 6 and 20 minutes for Western Blot analysis w/ anti-GST + anti-FLAG. in 1 ladder
- Western blots were each loaded w/ 5- μ L of each aliquot (4 μ L of protease soln. with 1 μ L of 5x SDS protein loading buffer.)

(see following page for anti-FLAG and anti-GST gels)

suggested: 1 μ g Factor Xa for 50 μ g protein

2 a,	$t=0$
3 a,	$t=1$ hr
4 a,	$t=3$ hr
5 a,	$t=6$ hr
6 a,	$t=20$ hr
7 b,	$t=0$
8 b,	$t=1$ hr
9 b,	$t=3$ hr
10 b,	$t=6$ hr
11 b,	$t=20$ hr
12 c,	$t=0$
13 c,	$t=1$ hr
14 c,	$t=6$ hr
15 c,	$t=20$ hr

Factor Xa cleavage at 4°C.

- 118 μL rxns set up (before Fact Xa added), then 24 μL removed (added to 24 μL 2x SDS loading buffer) and Fact. Xa added:


12 μL Factor Xa buffer
 48 μL (24 μg) pure protein (fr. p. 39)
 58 μL DI water
 (24 μL of above soln. removed, then)
 1 μL Fact. Xa dilution (0.1 and 0.013 μg)

- The rxn was incubated at 4°C. 24- μL aliquots were removed after 5 hours and after 30 hours. (Aliquots boiled w/ 2x SDS loading buffer. 20 μL of mixture loaded for Coomassie and 8 μL for the Western.)

anti-FLAG




1 2 3 4 5 6 7



7 6 5 4 3 2 1

in 1 ladder
 2 \rightarrow 0.1 μg per 20 μg protein, $t=0$
 3 " " $t=5$ hr
 4 " " $t=30$ hr
 5 \rightarrow 0.013 μg Fact Xa per 20 μg protein, $t=0$
 6 " " $t=5$ hr
 7 " " $t=30$ hr

Factor Xa cleavage with the addition of ZnCl




10 9 8 7 6 5 4 3 2 1
 (5 μ L protein soln. and 5 μ L of
 2 \times SDS loading buffer per lane.)

4 μ M GST + FLAG pure protein
 600 μ M ZnCl

20 μ L (20 μ g) protein
 3 μ L ZnCl soln. (3 μ g, 150 eq.)
 5 μ L 10 \times Factor Xa buffer
 31 μ L water
 (10 μ L of above soln. removed for
 t=0 point before Factor Xa added)
 1 μ L Factor Xa (0.5 or 0.25 ~~to~~ units)

lanes 2-4 0.5 units Novagen Factor Xa
 lanes 5-10 0.25 units Novagen Factor Xa
 (rxn. run at 24 $^{\circ}$ C)



1 2 3 4 5 6 7
 (same time course as above gel)

4 μ M pure protein
 200 μ M ZnCl

20 μ L (20 μ g) protein
 1 μ L ZnCl (1 μ g, 50 eq.)
 5 μ L Factor Xa 10 \times buffer
 33 μ L water
 (10 μ L of above soln. removed for t=0)
 1 μ L Factor Xa (0.5 or 0.25 u)

lanes 2-4 0.5 units Factor Xa
 lanes 5-7 0.25 units Factor Xa
 (rxn run at 16 $^{\circ}$ C)

Grr! The hope was that adding ZnCl might improve folding and thus decrease secondary site cleavage by organizing the four LIM domains (double Zn fingers - nm binding - that bind a total of 8 Zn.)

2 $^{\circ}$ cutting seems to occur faster than at the 1 $^{\circ}$ site to give fragments: \approx 43 kDa, 30 and 34 kDa and $<$ 10 kDa.
 There are no other IEG-R (Factor