



Tuesday, October 18, 2004

- The 16-mL GST-purified protein solution in TBS was incubated with 1.4 mL of FLAG-affinity resin (used 1x previously) w/ 40 mL of protease inhibitor cocktail at rt for 1.5 hours (then stored in fridge for 1 hour while at a meeting.)
- The resin was filtered, washed with 80 mL TBS, and then eluted with 1 to 1.2 mL aliquots of 0.1M glycine, pH 3.5 collected into 40 mL of 500 mM Tris·HCl, pH 8.0. The concentration of protein/sample was determined using the Biorad protein assay. All fractions were stored at 4°C.

Wednesday, Oct. 19, 2004

TEV cleavage (2 rxns run currently)

- Fraction 3 (600 µg in 1.1 mL) was split into 2 portions of 300 µg GST-ENLYFQC-Pax (38-557)-FLAG in 550 µL elution buffer/Tris. The following rxn was incubated at 37°C for 4 hours:
  - ≈ 300 µg protein in 550 µL
  - 60 µL of 500 mM Tris·HCl, pH 8.0
  - 15 µL mTEVp



1° 1.5 mL/10 mL TBST anti-FLAG ab

- | Lane | Description                              |
|------|--|
| 1    | MW marker                                |
| 2    | FLAG elution 1 (↓ Biorad quantification) |
| 3    | FLAG el. 2 (60 µg/mL)                    |
| 4    | FLAG el. 3 (600 µg/mL)                   |
| 5    | FLAG el. 4 (400 µg/mL; 1.1 mL)           |
| 6    | FLAG el. 5 (180 µg/mL; 1.1 mL)           |
| 7    | FLAG el. 6                               |
| 8    | TEV cleavage A, t = 4 h                  |
| 9    | TEV cleavage B, t = 4 h                  |

- The mTEVp cleavage reactions were stored at 4°C overnight, then used in ligation reactions.

EPL: pY31

Wednesday, Oct. 19, 2004




- 200  $\mu$ g of Cys(Pax38-557)-FLAG (from TEVp cleavage A on p. 65) in  $\approx$  600  $\mu$ L of mTEVp cleavage mixture was concentrated to 60  $\mu$ L using a 50-kDa MWCO centrifugal filter.

- 10  $\mu$ L were removed for gel analysis. The remaining 50  $\mu$ L were added to 250  $\mu$ g of lyophilized Achis6-Pax(2-36, pY31)-SBr and aliquots were removed and reacted as follows:

AII: 20  $\mu$ L of protein/thioester, 10  $\mu$ L of 50 mM Tris, 300 mM NaCl, pH 8.0, 2  $\mu$ L of 1 M MESNA. (No precip. after 24 h)

AIII: 10  $\mu$ L protein/thioester soln., 10  $\mu$ L of 200 mM MES, 300 mM NaCl, pH 6.1, 2  $\mu$ L of 1 M MESNA (100 mM final concentration) (A tiny amount of precipitation after 24 h.)

AIII: 10  $\mu$ L protein/thioester, 10  $\mu$ L of 200 mM MES, 300 mM NaCl, pH 6.1, 0.8  $\mu$ L of 1/20 dil. of thiophenol (0.2%) (some precip. after 24 h.)



1 2 3 4 5 6 7

5 4 3 2 1

5 4 3 2 1

\*

\*

\*

1° ab anti-His6

1° ab rabbit anti pY31-Pax

(5  $\mu$ L in 10 mL TBST)

(1.5  $\mu$ L in 10 mL TBST)

In 1 benchmark ladder

2 10  $\mu$ L Cys-Pax(38-557)-FLAG

3 rxn A, t=24 h

4 rxn AII, t=24 h, supernatant only (a tiny amnt of pellet got taken up also)

5 " , t=24 h, 1/2 of pellet resuspended

6 rxn AIII, t=24 h, supernatant only

7 " , t=24 h, 1/2 of pellet resuspended

(In each reaction  $\approx$  2.5 mg/mL thioester. - 0.5 mM)