

Phosphorylation of GST-Pax(38-557)-FLAG and Ac-His \times 6-Pax(2-557)-FLAG

Wednesday, July 27, 2005

Substrates used:

- GTP: Elution 4 fr. p. 137 of GST-Pax(38-557)-FLAG : 265 μ g/mL, 3.77 μ L/ μ g
 LigP: 1/3 dil. of NCL product from p. 136 (assuming 1.5 mg/mL original concentration from an ideal [I] of 3): 500 μ g/mL, 2 μ L/ μ g
 ACH \times 6-Pax(2-557)-FLAG

- Src, ERK, and JNK2 (Biosource phosphorylation) $\uparrow\uparrow$ JNK2 concentration
 3 40- μ L reactions were run w/ 2 μ g (7.6 μ L) GTP
 2 40- μ L reactions were run w/ 2 μ g (4 μ L) LigP:

2 μ g substrate (ie. 7.6 μ L)

8 μ L of 1 mM ATP (diluted from 10 mM in TBS) final [I] = 200 μ M

2.25 μ L of 226 mM MgCl₂ (diluted in TBS) final [I] = 15 mM

1.5 μ L Src (450 units), or 4 μ L ERK (148 units), or 7 μ L JNK2

TBS to 40 μ L total (20.65, 18.15, and 15.15 for Src, ERK, + JNK2)

The above reagents were combined in a 0.65 mL eppendorf, pipetted 7 \times to mix, flicked once and incubated at 30 $^{\circ}$ C for 15 min. 40 μ L of 2 \times SDS reducing gel loading buffer was added for a final substrate concentration of 25 ng/ μ L. A 1/5 dil. (10 μ L in 40 μ L 1 \times loading buffer) was made to give a 5 ng/ μ L solution.

- JNK1 (upstate) phosphorylation of GST-Pax(38-557)-FLAG (GTP)

2 50- μ L reactions:

2 μ g substrate (7.6 μ L)

8 μ L of 1 mM ATP

2.25 μ L of 226 mM MgCl₂

4 μ L of 0.1% Brij

4 μ L of 10 mM dtt

4 μ L of 50 mM EGTA

TBS to 50 μ L total (10.15 or 15.15 μ L)

5 μ L (0.5 μ g) or 10 μ L (1 μ g) of JNK1 α 1

The mixtures were pipetted 7 times, flicked, then incubated at 30 $^{\circ}$ C for 15 min. 30 μ L of 2 \times loading buffer were added (to give 25 ng/ μ L) and a 1/5 dilution was made (to give 5 μ L/ μ L).

Chemiluminescent gels from p. 138 phosphorylation (pan, pY31, and pY118)

Thursday, July 29, 2005


1° mouse anti-paxillin (1/1,000)

2° goat anti-mouse HRP (1/10,000)

exposed for 4 min.

GTP = GST-Pax(38-557)-FLAG

LigP = AcHisX6-Pax(2-557)-FLAG



- In
- 1 7 mL ladder
 - 2 JNK1 α 1 (0.5 mg Kinase in rxn)
 - 3 JNK1 α 1 (1.0 mg kinase in rxn)
 - 4 JNK2 (7 μ L!) for GTP
 - 5 GTP - no kinase
 - 6 LigP - no kinase
 - 7 ERK for GTP Src
 - 8 ERK for GTP
 - 9 Src for LigP
 - 10 ERK for LigP


30 ng loaded for anti-pax, pY31 and pY118
 250 ng loaded for anti-pS128, pS178

6 μ L (of 1/5 dilution) - 30 ng loaded/well

1° rabbit anti-pY31 pax (1/2,000)

2° goat anti-rabbit HRP (1/20,000)


exposed for 5 min.



1° rabbit anti-pY118 pax (1/2,000)

2° goat anti-rabbit HRP (1/20,000)

exposed for 3 minutes

30 ng (6 μ L of a 1/5 dilution) loaded per well

Chemiluminescent detection of pY31 Native Chem. Lig. product

Monday, August 29

- 1° rabbit anti-pY31 pax (1,2,000)
 2° goat anti-rabbit HRP (1,20,000)




- In 1 7 mL protein benchmark
 2 3 mL pY31 B (1/5 of A)
 3 6 mL pY31 B
 4 12 mL pY31 B
 5 24 mL pY31 B
 6 5 mL pY31 soln. A
 7 10 mL pY31 A
 8 15 mL pY31 A
 9 20 mL pY31 A

It appears that loading between 5 and 10 mL of the pY31 sample B results in the best visualization.

9 8 7 6 5 4 3 2 1
 expose time = 30 seconds
 run on a 0.75 width gel

Stock gel solutions for pY31 paxillin

130 mL (of the 300 mL total on p.146) with 130 mL 2x SDS page



- In 1 ladder
 2 30 mL FLAG (stock soln.)
 3 before ligation (8 mL + 22/30)
 4 15 mL FLAG
 5 NCL w/ pY31 (30 mL)
 6 NCL w/ pY31 (30 15 mL)
 7 NCL Y31 (7/20)
 8 NCL glycerol stock