


11-7-76

actually 11-8-76

3. Preparation of Crystalline Adduct

Sample: 11-3-1 Sample had a copious ppt  
cent down in several runs on PR-6 at  
2000 RPM in 30 ml Corex tubes.  
- Initial vol = 76 ml

4. Preparation of More Adduct Crystals

This sample will be used in the acetylation series to be run tomorrow.

- Crystals were taken from 10-26-2 SN which is the crystals from the mother liquor of the 1<sup>st</sup> large batch of adduct.
- These crystals are 11-7-5-J